



Injection Vials ISO 8362-1

Type		Neck Finish
2R	35 x 16 mm 1,0 mm	13 mm
4R	45 x 16 mm 1,0 mm	13 mm
6R	40 x 22 mm 1,0 mm	20 mm
8R	45 x 22 mm 1,0 mm	20 mm
10R	45 x 24 mm 1,0 mm	20 mm
15R	60 x 24 mm 1,0 mm	20 mm
20R	55 x 30 mm 1,2 mm	20 mm
25R	65 x 30 mm 1,2 mm	20 mm
30R	75 x 30 mm 1,2 mm	20 mm

other sizes on request
glass type 1 clear or amber
glass type 2 amber
glass type 3 clear