

Large Volume Vials

MG Standard Infusion Vials



Volume	Dimension	Neck Finish
75 ml	75,0 x 47,0 / 1,95	28 mm
100 ml	95,0 x 47,0 / 1,95	28 mm
125 ml	115,0 x 47,0 / 1,95	28 mm
75 ml	75,0 x 47,0 / 1,95	32 mm
100 ml	95,0 x 47,0 / 1,95	32 mm
125 ml	115,0 x 47,0 / 1,95	32 mm

other sizes on request
glass type 1 clear
glass type 1 amber on request

Quantity per Tray
Vial Diameter 47,0 mm 67 pcs.

Tray Dimension: 392 x 370 x 60 mm
Minimum Order Quantity: 3000 pcs., exw



Münnerstädter Glaswarenfabrik GmbH
Otto-Liebmann-Straße 2
DE-97702 Münnerstadt
Phone: +49 (0) 97 33 / 8 08-0
Fax: +49 (0) 97 33 / 8 08-160
Internet: <http://www.mglas.de>
E-Mail: info@mglas.de

