



## Large Volume Vials

### MG Standard Injection Vials

| Volume | Dimension          | Neck Finish |
|--------|--------------------|-------------|
| 50 ml  | 73,0 x 40,0 / 1,5  | 20 mm       |
| 75 ml  | 100,0 x 40,0 / 1,5 | 20 mm       |
| 100 ml | 120,0 x 40,0 / 1,5 | 20 mm       |
| 75 ml  | 75,0 x 47,0 / 1,7  | 20 mm       |
| 100 ml | 95,0 x 47,0 / 1,7  | 20 mm       |
| 125 ml | 115,0 x 47,0 / 1,7 | 20 mm       |

**other sizes on request**  
**glass type 1 clear**  
**glass type 1 amber on request**

**Quantity per Tray**  
**Vial Diameter 40,0 mm 90 pcs.**  
**Vial Diameter 47,0 mm 67 pcs.**

**Tray Dimension: 392 x 370 x 60 mm**  
**Minimum Order Quality: 3000 pcs., exw**



Münnerstädter Glaswarenfabrik GmbH  
Otto-Liebmann-Straße 2  
DE-97702 Münnerstadt  
Phone: +49 (0) 97 33 / 8 08-0  
Fax: +49 (0) 97 33 / 8 08-160  
Internet: <http://www.mglas.de>  
E-Mail: [info@mglas.de](mailto:info@mglas.de)

