



Large Volume Vials

MG Standard Screw Neck Vials

| Volume | Dimension | Neck Finish |
|--------|--------------------|---------------|
| 50 ml | 70,0 x 40,0 / 1,5 | 22 |
| 60 ml | 80,0 x 40,0 / 1,5 | 22, 28, 32 |
| 70 ml | 90,0 x 40,0 / 1,5 | 22, 28, 32 |
| 75 ml | 95,0 x 40,0 / 1,5 | 22, 28, 32 |
| 80 ml | 100,0 x 40,0 / 1,5 | 22, 28, 32 |
| 90 ml | 110,0 x 40,0 / 1,5 | 22, 28, 32 |
| 100 ml | 120,0 x 40,0 / 1,5 | 22, 28, 32 |
| 75 ml | 80,0 x 47,0 / 1,7 | 22, 28, 32,40 |
| 80 ml | 85,0 x 47,0 / 1,7 | 22, 28, 32,40 |
| 90 ml | 90,0 x 47,0 / 1,7 | 22, 28, 32,40 |
| 100 ml | 95,0 x 47,0 / 1,7 | 22, 28, 32,40 |
| 125 ml | 110,0 x 47,0 / 1,7 | 22, 28, 32,40 |

other sizes on request
glass type 1 clear
glass type 1 amber on request

Quantity per Tray
Vial Diameter 40,0 mm 90 pcs.
Vial Diameter 47,0 mm 67 pcs.

Tray Dimension: 392 x 370 x 60 mm
Minimum Order Quality: 3000 pcs., exw



Münnerstädter Glaswarenfabrik GmbH
Otto-Liebmann-Straße 2
DE-97702 Münnerstadt
Phone: +49 (0) 97 33 / 8 08-0
Fax: +49 (0) 97 33 / 8 08-160
Internet: <http://www.mglas.de>
E-Mail: info@mglas.de

