



Screw Neck Vials ISO 11418-7

Capacity	Nominal Thread
5 ml 45 x 18 mm 1,2	14 x 2,5
7,5 ml 40 x 22 mm 1,2	18 x 3
10 ml 45 x 24 mm 1,2	18 x 3
15 ml 60 x 24 mm 1,2	18 x 3
20 ml 55 x 30 mm 1,2	22 x 3
25 ml 65 x 30 mm 1,2	22 x 3
30 ml 75 x 30 mm 1,2	22 x 3

other sizes on request
glass type 1 clear or amber
glass type 2 amber
glass type 3 clear