












Pre-fillable Glass Syringes ISO 11040-4






other sizes on request
glass type 1 clear

Münnerstädter Glaswarenfabrik GmbH
Otto-Liebmann-Straße 2
DE-97702 Münnerstadt
Phone: +49 (0) 97 33 / 8 08-0
Fax: +49 (0) 97 33 / 8 08-160
Internet: <http://www.mglas.de>
E-Mail: info@mglas.de



	Needle Shield (NS)	Needle (N)	Glassbarrel (B)	Plunger Stoppers (PS)	Plunger Rods (PR)
					
Syringe Type	NS	N	B	PS	PR
0,5 ml 27 G 1/2	3003	0,40 x 12,7	Ø 6,85 x 57,5	2125	E258-01.4
1,0 ml 27 G 1/2 long	3003	0,40 x 12,7	Ø 8,15 x 64,0	2243	E216-02.4
1,0 ml 27 G 1/2	3003	0,40 x 12,7	Ø 10,85 x 46,7	2116	E142-01.4
1,0 ml 25 G 5/8	3005	0,35 x 16,0	Ø 10,85 x 46,7	2116	E142-01.4
1,0 ml 23 G 1	3009	0,60 x 25,5	Ø 10,85 x 46,7	2116	E142-01.4

	Tip Cap (TC)	Glassbarrel (B)	Plunger Stoppers (PS)	Plunger Rods (PR)
				
Syringe Type	TC	B	PS	PR
0,5 ml Luer	3379	Ø 6,85 x 58,6	2125	E258-01.4
1,0 ml Luer long	3379	Ø 8,15 x 64,0	2243	E216-02.4
1,0 ml Luer	3379	Ø 10,85 x 46,7	2116	E142-01.4
2,25 ml Luer	3379	Ø 10,85 x 65,5	2116	E146-01.4
3,0 ml Luer	3379	Ø 10,85 x 83,2	2116	E147-01.4
5,0 ml Luer	3379	Ø 14,45 x 79,4	2114	E148-01.4

	Tip Cap (TC)	MG Luer-Lock Adapter (LAL)	Glassbarrel (B)	Plunger Stoppers (PS)	Plunger Rods (PR)
					
Syringe Type	TC	LAL	B	PS	PR
0,5 ml Luer-Lock	3379	839127	Ø 6,85 x 61,1	2125	E258-01.4
1,0 ml Luer-Lock long	3379	839127	Ø 8,15 x 65,0	2243	E216-02.4
1,0 ml Luer Lock	3379	839127	Ø 10,85 x 48,3	2116	E142-01.4
2,25 ml Luer-Lock	3379	839127	Ø 10,85 x 67,0	2116	E146-01.4
3,0 ml Luer-Lock	3379	839127	Ø 10,85 x 85,7	2116	E147-01.4
5,0 ml Luer-Lock	3379	839127	Ø 14,45 x 80,0	2114	E148-01.4